
Instructions for use

Miniflow

Contents

EN

DE

CZ

DA

EL

ES

FI

FR

HU

IT

JA

NL

NO

PL

PT

RO

RU

SL

SV

TR

The user has to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State/Country in which the user is established.

Serious incidents are defined as malfunctions, failures or changes in characteristics or performance, or inadequacies in the labelling or instructions for use of a medical device, which directly or indirectly have led, could have led or could lead to the death or serious deterioration of the state of health of a patient, user or other person.

For Masks:

CE mark according to the current conformity assessment procedure and declaration of conformity (see website).

Classification (MDR (EU) 2017/745 Annex VIII): I

For Miniflow and Prongs:

CE mark according to the current conformity assessment procedure and declaration of conformity (see website).

Classification (MDR (EU) 2017/745 Annex VIII): Iia

medin Medical Innovations GmbH
Adam-Geisler-Str. 1 - 82140 Olching - Germany
+49 8142 448 460 - www.medin-medical.com - info@medin-medical.com
Made in Germany

IFU_Miniflow_Rev12_EN

Miniflow:

REF 4000, 4001, 6684635, BM4000SLE, 10162671, 10162673, 10162674, 4000-10, 6684897, 6684900, BM4000/10SLE

Prongs:

REF 1200-01, 1200-02, 1200-03, 1200-21, 1200-22, 1200-32, 1200-33, 10162662, 6684645, 6684647, 6684650, 6684652, 6684655, 6684657, 6696610, 10162664, BM1026, BM1200/02, BM1200/03, BM1200/21, BM1200/22, BM1200/32, BM1200/33, 10162665

Masks:

REF 1200-04, 1200-05, 1200-06, 1200-07, 1200-08, 1200-09, 10162666, 6684637, 6684640, 6684642, 6688545, 10162668, BM1200/04, BM1200/05, BM1200/06, BM1200/07, BM1200/08, 10162669

© 2023 medin Medical Innovations GmbH. All rights reserved. Printed in Germany.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of medin Medical Innovations GmbH.

This document may be revised, replaced, or made obsolete by other documents by medin Medical Innovations GmbH at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the customer support department of medin Medical Innovations GmbH or download the current applicable version from our homepage www.medin-medical.com. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way the right of medin Medical Innovations GmbH to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, medin Medical Innovations GmbH has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. medin Medical Innovations GmbH sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in the user manual.

medin Medical Innovations GmbH shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-medin medical Innovations GmbH parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

medin Medical Innovations GmbH will make available, on request, circuit diagrams, component parts lists, descriptions, calibration instructions, or other information that will assist appropriately trained personnel to repair those parts of the equipment designated by medin Medical Innovations GmbH to be repairable.

For all proprietary as well as third-party trademarks used by medin Medical Innovations GmbH, see www.medin-medical.com/de/Patents.html. Product and/or company names marked with a may be the trademarks and/or registered trademarks of their respective owners.

Table of contents

1	Introduction.....	
1.1	About these instructions for use	
1.2	General liability conditions	
2	Intended use.....	
2.1	Miniflow.....	
2.2	Prong	
2.3	Mask	
2.4	Indications, expected clinical benefit, contraindications	
3	Principles of function and technical specifications.....	
3.1	Miniflow.....	
3.2	Bonnets, prongs, masks, Minifoam	
3.2.1	Bonnets.....	
3.2.2	Prongs	
3.2.3	Masks.....	
3.2.4	Minifoam	
3.3	Environmental conditions (operation / transport / storage).....	
4	Assembly and use	
4.1	Equipment needed	
4.2	Additional aids	
4.3	Structure and set-up of the system	
5	Products which can be connected to the Miniflow.....	
5.1	Miniflow.....	
5.2	Prongs, masks	
6	Accessories	
7	Disposal	
8	Symbols	
9	Abbreviations	

1 Introduction

1.1 About these instructions for use

These instructions for use contain information on starting up and using the Miniflow/prong/mask. They also contain safety information, describe the product function, and give an overview of the accessories needed.

Users must be very familiar with the information and warnings given in these instructions for use in order to use the Miniflow/prong/mask safely. Nonetheless, they are no substitute for training. The instructions and warnings are categorized as follows:

Warning: Warnings must be followed in order to prevent possible serious consequences for the patient or the user.

Attention: This indicates hazards which could damage the device or impair its functionality.

Note: This indicates how the Miniflow/prong/mask can be used more efficiently.

Keep these instructions for use close to the Miniflow/prong/mask and easily accessible.

Please contact the dealer or go directly to the manufacturer if you have any questions or comments about these instructions for use.

Warning:

- Read these instructions for use carefully before using the products, and follow the instructions, warnings and tips closely.
- The Miniflow/prong/mask may only be used by individuals with a thorough knowledge of and in accordance with these instructions for use for the purpose described under the section headed "Intended use".
- The use of devices and consumables used in combination with the Miniflow/prong/mask is described in the respective instructions for use. Users must also comply with these instructions for use.

1.2 General liability conditions

The general terms of trade of medin Medical Innovations GmbH (hereinafter referred to as medin) are binding.

medin accepts no responsibility or liability for the safe use of the Miniflow/prong/mask if an operation that does not comply with intended use is performed. The user is responsible for respecting the specifications in the instructions for use. In addition, the use of the Miniflow/prong/mask requires training of the user by a qualified person.

The duration of use of the Miniflow/prong/mask is defined as a maximum of 7 days.

medin accepts no liability for damage caused by failure to comply with these instructions.

2 Intended use

2.1 Miniflow

The miniflow is a single-use product for non-invasive respiratory support and oxygen therapy of newborn infants. It is used in combination with a ventilator or bubble CPAP device which supplies the necessary gas flow to generate the CPAP pressure. The miniflow is used under the supervision of competent personnel trained for this purpose in a clinical environment.

Warning:

- The oxygen saturation in the patient's blood must be continuously monitored during use.
- The Miniflow, masks, prongs, bonnets and Minifoam are disposable products; they may be used only on a single patient – and may not be reprocessed.
- The Miniflow may only be used in combination with prongs, masks, bonnets and Minifoam
- The duration of use of the Miniflow/prong/mask is defined as a maximum of 7 days.
- The product must be disposed of in accordance with the manufacturer's instructions

2.2 Prong

The prong is a single-use product for non-invasive respiratory support and oxygen therapy of newborn infants. It is used in combination with a medijet or miniflow. The prong is used under the supervision of competent personnel trained for this purpose in a clinical environment.

2.3 Mask

The mask is a single-use product for non - invasive respiratory support and oxygen therapy of newborn infants. It is used in combination with a medijet or miniflow. The mask is used under the supervision of competent personnel trained for this purpose in a clinical environment.

2.4 Indications, expected clinical benefit, contraindications

Indications

Neonatal respiratory insufficiency

Contraindications

Non-invasive respiratory support is contraindicated if certain clinical situations, anatomic malformations, or the degrees of severity of a disease necessitates intubation and mechanical ventilation or make non-invasive treatment impossible

Benefit

The miniflow/prong/mask improves or maintains the gas exchange in the lungs of the patients

Intended patient / target group

Neonatal patients

Intended user

Miniflow/prong/mask is intended to be use by trained health care professionals under the direct supervision of a licensed physician

Intended environment

In a clinical environment

3 Principles of function and technical specifications

3.1 Miniflow

- The Miniflow is a patient interface which can be used in CPAP mode or during non-invasive ventilation in combination with a ventilator or a bubble CPAP system tested by medin for compatibility. The Miniflow can be connected via a double tubing system to a CPAP driver or ventilator. The CPAP parameters are adjusted and controlled on the ventilator. The ventilator also provides the necessary gas flow.
- In the Miniflow, the expiration and the inspiration gas flow are separated from each other by a partition until they reach the prong adapter. The prong adapter is tiltable and this allows the direction of the prong adapter to be adjusted to the baby's nasal anatomy and size which facilitates fixation.
- A prong or a mask selected according to the baby's size is placed on the prong adapter.
- The prong or the mask ensures a soft contact between the Miniflow and the baby's nose and can be used in alternation to avoid pressure points.
- The Miniflow itself is secured to the bonnet which the baby wears.

Technical data	
General	
Expiration connection	10 mm
Inspiration connection	10 mm
Weight	10 g
Duration of use	Max. 7 days
Parameters & sensors	
CPAP pressure	0 to > 20 mbar
Principles of function	Continuous flow
Connections	
Inspiration	OD: 9.9 – 10.9 mm, conical / ID: 7.5 mm
Expiration	OD: 9.9 – 10.9 mm, conical / ID: 7.5 mm

3.2 Bonnets, prongs, masks, Minifoam

3.2.1 Bonnets

A divided fleece is sewn onto the upper side of the bonnets. A foam wedge (Minifoam) is placed between the fleece parts in which the Miniflow is placed and secured with a separate Velcro strap. The straps additionally included are used to stabilize the mask or prong position and, with their Velcro ends, they can be placed anywhere on the bonnet fabric. The masks and prongs each contain loops for the fixation straps of the bonnet. In addition, they contain a pair of holes for the connection to the Miniflow in which the prong adapter of the Miniflow is inserted.

3.2.2 Prongs

The binasal prongs of medin (see Figure 1) serve as an interface between the Miniflow and the patient's nasal area and thus complete the nCPAP system. They are introduced into the nostrils and thus establish the connection to the system. The prongs are available in seven sizes: XS, S, M, MW, L, LW, XL. Two short tubes are attached to a base piece. On the side

of the base piece are fixation tabs through which the fixation straps of the bonnets are pulled in order to secure the prong to the bonnet. The base of the tubes has a tapered shape in order to seal the transition between the nostril and the prong through more flexible adjustment.


Figure 1 - Prong with fixation tabs (1) and tapered tubes (2)

The various prong sizes differ in diameter and distance between the tubes. The appropriate prong can be selected for the patient using a special measuring tape, in order to keep skin damage on the patient to a minimum. The tabs of the prong are attached at a slanted angle. This ensures a comfortable position on the patient's cheek during fixation to the bonnet and additionally prevents the tubes from spreading when the tabs are firmly pulled apart. The prongs can be used on the Miniflow and also on the Medijet.

3.2.3 Masks

Binasal masks (see Figure 2), like the binasal prongs, serve as an interface between the Miniflow and the patient's nose. medin offers masks in five sizes: XS, S, M, L and XL. They are placed on the patient's nose and are intended to form a tight seal. The masks have a malleable cushion which establishes contact with the skin. This cushion can compensate for uneven areas and because of its malleability, it reduces skin damage. The fixation straps are secured in the tabs attached to the side and the straps are then secured to the bonnet. The masks can be used on the Miniflow and also on the Medijet.


Figure 2 - Mask with connection opening (1) for the nCPAP interface, fixation tabs (2) and application cushion (3)

3.2.4 Minifoam

The Minifoam is a foam wedge which enables the Miniflow to be embedded and secured without additional aids. The Miniflow is fit into the Minifoam, see Figure 3. The opening in the Minifoam corresponds to the diameter of the Miniflow tube. The Minifoam is secured to the bonnet together with the Miniflow. The Velcro strap corresponds to the dimensions of the Velcro pads on the bonnet.


Figure 3 - Minifoam (1) and Miniflow fitted into Minifoam (2)

3.3 Environmental conditions (operation / transport / storage)

Environmental conditions for Miniflow, prong, mask during operation:

Temperature:	15°C to 40°C
Cleanliness:	The hygiene regulations valid for the hospital are to be observed. In addition, the Miniflow/prong/mask is to be used under standard hospital ambient conditions

Environmental conditions for Miniflow, prong, mask during storage/transport:

Temperature:	-30°C to 70°C
Cleanliness:	The Miniflow/prong/mask must be protected from contamination during transport and storage
Other:	The Miniflow/prong/mask must be stored dry and protected from sunlight

4 Assembly and use

4.1 Equipment needed

To perform effective CPAP therapy with the Miniflow, the following equipment is needed:

- A ventilator or a bubble CPAP system
- An active respiratory gas humidifier
- A Miniflow
- A Minifoam
- An appropriate double tubing system which must contain an inspiration and expiration line. In addition, a chamber for the humidifier is needed. The connections to the tubing ends must be F10 to fit on the Miniflow.
- A bonnet – its size must be selected such that it fits the baby.
- A mask or a prong which fits the size of the patient's nose
- Optional for ventilators which need external pressure measurement:
A T piece for measuring the pressure (already included in the medin tubing or can be ordered separately, REF 4010)

4.2 Additional aids

- A measuring tape to help select the sizes of bonnets, masks and prongs.

Warning:

If the packaging of a part is damaged or if the products themselves are damaged, these parts should not be used and should be disposed of.

If the openings are not patent, the products in question should not be used.

4.3 Structure and set-up of the system


Figure 4 - Set-up of the Miniflow, patient side


Figure 5 - tube system in combination with Miniflow

(1) Connect the tube set to the nCPAP driver or the ventilator and the humidifier and follow the installation instructions of the humidifier. If necessary, use the enclosed adapter.

(2) Connect the Miniflow to the inspiration and expiration tubing. If your ventilator requires external pressure measurement, insert the T piece for pressure measurement in the expiration part. The pressure measurement line can then be inserted in this piece and connected with the ventilator.

(3) Close the prong adapter of the Miniflow manually and adjust the CPAP parameters, such as the pressure, on the CPAP driver or ventilator. If this is not possible, check the tubing connections for leaks. If all of the connections are correct and the Miniflow is still causing problems, it should not be used and must be disposed of.

(4) Select a prong or a mask, using the measuring tape as an aid, if necessary (see Figure 6). The prong should be selected such that it is large enough to seal the entire nostril. The mask should be selected such that it fits precisely along the sides of the nostrils. A prong is placed on the Miniflow such that the arched side of the prong points toward the patient's upper lip. The mask is attached such that it conforms to the shape of the nose.


Figure 6 - Determining the size of the head (a) or the prong/mask size (b) using the measuring tape

Note:

The mask size should be selected such that it fits precisely along the sides of the nostrils. The mask is attached such that it conforms to the shape of the nose.

Warning:

In addition, check the holes of the Miniflow and the prongs/masks to ensure they are patent.

Note:

Always use the largest possible prong to avoid leaks and ensure the best possible positioning – the smaller the leak is, the higher and more stable the CPAP pressure is.

Warning:

Always insert the Miniflow fully into the prongs or masks and do not use any cream near the prongs, masks or Miniflow since otherwise the prong or mask can slip from the Miniflow.

Attention:

The prongs and masks should not be reprocessed or reused. They are only for single use. Reprocessed prongs and masks may harden and injure the patient.

(5) Select a bonnet by means of the color coding and use the measuring tape as an aid, if necessary (see Figure 6). The bonnet should be secured firmly enough but without leading to deformation of the head. Remove the separate Velcro strap from the bonnet and put it aside, within easy reach. The bonnet is placed on the head such that the gap between the two fleece parts is centered over the forehead, without exerting too much pressure on the head. The bonnet should cover the ears and be pulled down to the nape of the neck (see Figure 7). Guide the fixation straps through the loops on the prong or on the mask.

Figure 7 - Positioning the bonnet on the patient's head

Note: Bonnets

Pull the bonnet down to just above the patient's eyebrows, as close

as possible to the eyes – this makes it easier to secure the Miniflow. Bonnets which are too small lead to deformation of the patient's head. Bonnets which are too large cannot ensure sufficient fixation of the Miniflow.

(6) Select the most appropriate angle for the prong adapter of the Miniflow (45° or 60°), see Figure 8. If a prong is used, the narrower angle is generally more suitable; if a mask is used, the wider angle. Place the Minifoam foam wedge between the two fleece pieces in the middle of the bonnet (see Figure 9). Fit the Miniflow into the oval opening of the foam provided. Insert the prong into the nostrils. In doing so, maintain a small distance between the nasal septum and the base of the prong. When a mask is used, ensure that it fits evenly on the skin and that a distance is maintained between the inner rim of the mask and the nostrils. Using the Velcro strap, secure the Miniflow from above to the fleece pieces of the bonnet. Secure the Velcro area of the straps in the area of the ears as an extension of the fixation tabs of the prong. Do not pull the straps too tightly; they should only provide the lateral support needed.

Figure 8 - Adjusting the angle of the Miniflow

Figure 9 - Securing the Miniflow to the bonnet

Attention:

Ensure that you mount mask/prong in the correct way, see figure 8

Attention:

- Do not insert the prongs fully into the nostrils. The connector block should not touch the nose.
- Regularly switching between prong and mask is an effective method for avoiding or reducing damage to the patient's skin.
- The Miniflow is held primarily by means of the bonnet and the foam. If the leak is too large, check the selected prong or mask size – pulling the straps more tightly does not by itself solve the problem.
- Do not pull the straps too tight. Excessively tightening the straps or inserting the prongs too far can lead to pressure points on the baby's nose or skin.

Attention:

- Ensure that the weight of the tubes does not pull directly on the patient's head or on the interface. Use cushioning to support the tubes.

(7) Regularly check the position and fixation of the Miniflow and the accessories.

Attention:

Regularly check the positioning and fixation of the Miniflow and the accessories since pressure points may otherwise develop. Check the prongs and masks regularly for collected secretions since there otherwise may be an obstruction.

5 Products which can be connected to the Miniflow

5.1 Miniflow

The Miniflow shall only be combined with the following devices tested by medin for compatibility:

Product name	REF	Description
medinBLENDER	1090 1091 1093 1094	To provide the volume flow rate
Bubble Water Column with overpressure valve	2040	
10x nCPAP circuit set (dual-limb)	206746	For connection with the driver device
F&P 900MR850	MR850ARU	Humidifier Fisher&Paykel

The use of products that are not compatible puts the patient at risk and may impair the proper and correct functioning of the product. Combining the Miniflow with another product not listed falls under the responsibility of the user.

5.2 Prongs, masks

Prongs and masks shall only be combined with the Medijet or the Miniflow.

6 Accessories

medin provides various components as accessories for the use of the Miniflow:

Product name	REF	Comments
10x Prong (x-small)	1200-01	
10x Prong (medium)	1200-02	
10x Prong (x-large)	1200-03	
10x Prong (small)	1200-21	
10x Prong (large)	1200-22	
10x Prong (medium-wide)	1200-32	
10x Prong (large-wide)	1200-33	
Medin Mixed Box Prong (1XS, 2S, 2M, 1MW, 2L, 1LW, 1XL)	10162662	
10x Mask (small)	1200-04	
10x Mask (medium)	1200-05	
10x Mask (large)	1200-06	
5x Mask (x-large)	1200-07	
10x Mask (x-small)	1200-08	
10x Mask (large +)	1200-09	
Medin Mixed Mask (1XS, 1S, 1M, 1L, 1XL)	10162666	
10x Bonnet (xx-small)	1213-10	light green
10x Bonnet (x-small)	1214-10	white
10x Bonnet (small)	1215-10	yellow
10x Bonnet (medium)	1216-10	red
10x Bonnet (large)	1217-10	light blue
10x Bonnet (x-large)	1218-10	orange
10x Bonnet (xx-large)	1219-10	light green
10x Bonnet (xxx-large)	1220-10	white
10x Minifoam	4030	Fixation cushion for Miniflow

Warning: Accessories

The Miniflow only functions in combination with these masks, prongs, bonnets and the Minifoam, otherwise proper function cannot be guaranteed.

7 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

8 Symbols

The symbols explained below are used as labels on the product or in these instructions for use:

Symbol	Meaning	Symbol	Meaning
	CE mark		Temperature limit
	Medical device		Latex-free
	Follow instructions for use		Manufacturer
	Article number		Protect from sunlight
	Production lot number, batch		Date of manufacture
	Do not reuse		Use before
	Keep dry		Quantity
	For use by medical prescription only		Do not use if packaging is Damaged
	Authorized representative		Revision index
	Distributed by		

9 Abbreviations

EN

Abbreviation	Meaning
%	Percent
(m)bar	(Milli-)bar
(n)CPAP	(nasal) Continuous Positive Airway Pressure
°C	Degrees Celsius
g	Gram
ID	Inner diameter
mm	Millimeter
OD	Outer diameter

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© 2023 medin Medical Innovations GmbH. All rights reserved. ドイツ印刷