
Instructions for use

Nuflow

Contents

EN

DE

CZ

DA

EL

ES

FI

FR

HU

IT

JA

NL

NO

PL

PT

RO

RU

SL

SV

TR

ZH

The user has to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State/Country in which the user is established.

Serious incidents are defined as malfunctions, failures or changes in characteristics or performance, or inadequacies in the labelling or instructions for use of a medical device, which directly or indirectly have led, could have led or could lead to the death or serious deterioration of the state of health of a patient, user or other person.

CE 0483

For Nuflow

CE mark according to the current conformity assessment procedure and declaration of conformity (see website).

Classification (MDR (EU) 2017/745 Annex VIII): IIa

medin Medical Innovations GmbH

Adam-Geisler-Str. 1 - 82140 Olching - Germany

+49 8142 448 460 - www.medin-medical.com - info@medin-medical.com

Made in Germany

IFU_Nuflow_Rev05_EN

REF 1320, 1330, 1340, 1350, 1393, 1394, 10072354, 10072355, 10072356, 10072357

© 2023 medin Medical Innovations GmbH. All rights reserved. Printed in Germany.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of medin Medical Innovations GmbH.

This document may be revised, replaced, or made obsolete by other documents by medin Medical Innovations GmbH at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the customer support department of medin Medical Innovations GmbH or download the current applicable version from our homepage www.medin-medical.com. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way the right of medin Medical Innovations GmbH to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, medin Medical Innovations GmbH has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. medin Medical Innovations GmbH sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in the user manual.

medin Medical Innovations GmbH shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-medin medical Innovations GmbH parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

medin Medical Innovations GmbH will make available, on request, circuit diagrams, component parts lists, descriptions, calibration instructions, or other information that will assist appropriately trained personnel to repair those parts of the equipment designated by medin Medical Innovations GmbH to be repairable.

For all proprietary as well as third-party trademarks used by medin Medical Innovations GmbH, see www.medin-medical.com/de/Patents.html. Product and/or company names marked with a may be the trademarks and/or registered trademarks of their respective owners.

Table of contents

1	Introduction.....	
1.1	About these instructions for use	
1.2	General liability conditions	
2	Intended use.....	
2.1	Nuflow	
2.2	Indications, expected clinical benefit, contraindications	
3	Principles of function and technical specification	
3.1	Nuflow nasal cannula	
3.2	Environmental conditions (operation / transport / storage).....	
4	Assembly and use	
4.1	Equipment needed	
4.2	Structure and set-up of the system	
5	Products which can be connected to the Nuflow	
6	Accessories	
7	Disposal	
8	Symbols	
9	Abbreviations	

1 Introduction

1.1 About these instructions for use

These instructions for use contain information on starting up and using the Nuflow. They also contain safety information, describe the product function, and give an overview of the accessories needed.

Users must be very familiar with the information and warnings given in these instructions for use in order to use the Nuflow safely. Nonetheless, they are no substitute for training. The instructions and warnings are categorized as follows:

Warning: Warnings must be followed in order to prevent possible serious consequences for the patient or the user.

Attention: This indicates hazards which could damage the device or impair its functionality.

Note: This indicates how the Nuflow can be used more efficiently.

Keep these instructions for use close to the nCPAP circuit sets and easily accessible.

Please contact the dealer or go directly to the manufacturer if you have any questions or comments about these instructions for use.

Warning:

- Read these instructions for use carefully before using the products, and follow the instructions, warnings and tips closely.
- The Nuflow may only be used by individuals with a thorough knowledge of and in accordance with these instructions for use for the purpose described under the section headed "Intended use".
- The use of devices and consumables used in combination with the Nuflow is described in the respective instructions for use. Users must also comply with these instructions for use.

1.2 General liability conditions

The general terms of trade of medin Medical Innovations GmbH (hereinafter referred to as medin) are binding.

medin accepts no responsibility or liability for the safe use of the Nuflow if an operation that does not comply with intended use is performed. The user is responsible for respecting the specifications in the instructions for use. In addition, the use of the Nuflow requires training of the user by a qualified person.

The duration of use of the Nuflow is defined as a maximum of 7 days.

medin accepts no liability for damage caused by failure to comply with these instructions.

2 Intended use

2.1 Nuflow

The Nuflow nasal cannula is a single - use product intended for nasal high - flow therapy or oxygen therapy for newborn infants up to young children.

Warning:

- The Nuflow nasal cannula may be used only on a single patient and may not be reprocessed.
- Overpressure protection (e.g. an overpressure valve) must always be present in the system.
- The duration of use of the Nuflow nasal cannula is defined as a maximum of 7 days
- The Nuflow nasal cannula may only be used with the accessories specified in chapter 6
- The Nuflow nasal cannula may only be used in combination with the products tested for combinability specified in chapter 5
- The product may be disposed of in accordance with the manufacturer's instructions, see chapter 7 - Disposal.

2.2 Indications, expected clinical benefit, contraindications

Indications

Neonatal respiratory insufficiency.

Contraindications

High-flow therapy is contraindicated if certain clinical situations, anatomic malformations, or the degrees of severity of a disease necessitate intubation and mechanical ventilation or make high-flow treatment impossible.

Benefits

Nasal high-flow therapy improves or maintains the gas exchange in the lungs of the patients.

Intended patient / target group

Newborn infants up to young children

Intended user

The Nuflow nasal cannula is intended to be used by trained healthcare professionals under the direct supervision of a licensed physician.

Intended environment

In a clinical environment while monitoring the patient's oxygen saturation.

3 Principles of function and technical specification

3.1 Nuflow nasal cannula

The Nuflow nasal cannula is a patient interface (interface between the patient and a suitable high-flow device) for nasal high-flow therapy. To achieve this, the Nuflow nasal cannula consists of a pair of tubes which can be connected on one end with the ventilation tube and which end in a prong adapter on the other end which is positioned in the patient's nostrils.

The required flow and required oxygen concentration are adjusted on the high-flow device. If the Nuflow nasal cannula is used on a ventilator with High Flow mode, the adapter for the Y-piece (OD 15 mm/ID 10 mm) of the tubing system may be needed. The size of the Nuflow nasal cannula is selected according to the size of the patient's nose. The Nuflow nasal cannula is secured to the patient's skin either with the fixation patch provided or it must be secured directly to the skin using a suitable self-adhesive medium (e.g. fixation gauze/viscose gauze).

After the entire system has been set up (insertion of the Nuflow nasal cannula into the supplying respiratory tube), the application of the Nuflow nasal cannula in the patient's nose and its fixation to the patient's head, the Nuflow nasal cannula supplies the humidified and heated air/oxygen mixture transferred to it into the patient's nostrils.

A flushing effect is achieved there through the air flow into the nose and the dead space in the nose is reduced. Through the possible enrichment of inhaled air with oxygen (increase in FiO₂), the diffusion of the oxygen in the blood is increased and due to turbulence in the nostril, an increased pressure is generated in the nose.

Technical specifications	
Tube connections available	
OD 10 mm or OD 15 mm	
Weight of Nuflow nasal cannula incl. nasal bandage and adapter (without packaging)	
Nuflow nasal cannula size 1320 and 1330	13 g
Nuflow nasal cannula size 1340 and 1350	13 g
Sizes and flow limits of the Nuflow nasal cannula	
The following maximum flow (basis of measurement according to DIN 1343 based on 0°C, 1013.25 mbar and 0% air humidity (0/1013)) generates a backpressure of approx. 60 mbar in the tubing system. Whether this maximum flow can be achieved in practice also depends on the flow source used, the overpressure protection, leakage at the nostril, and the tubing system used.	
Nuflow nasal cannula (small)	Max. 10 L/min
Nuflow nasal cannula (medium)	Max. 14 L/min
Nuflow nasal cannula (large)	Max. 23 L/min
Nuflow nasal cannula (x-large)	Max. 27 L/min
Duration of use	
Max. 7 days	

3.2 Environmental conditions (operation / transport / storage)

Warning:

During transport, storage and operation, the environmental conditions specified in chapter 3.2 – Environmental conditions (operation /transport /storage) must be observed.

Environmental conditions for the Nuflow during operation:

Temperature:	0°C to 40°C
Cleanliness:	The hygiene regulations valid for the hospital are to be observed. In addition, the Nuflow is to be used under standard hospital ambient conditions.

Environmental conditions for the Nuflow during storage/transport:

Temperature:	0°C to 40°C
Cleanliness:	The Nuflow must be protected from contamination during transport and storage.
Other:	The Nuflow must be stored in a clean and dry room and protected from direct sunlight.

4 Assembly and use

4.1 Equipment needed

The high-flow therapy can be performed with various devices. The following equipment is needed in the configuration options below:

- **Combination with medinBLENDER/gas mixing unit**

- (1) Nuflow nasal cannula incl. fixation patch
- (2) medinBLENDER/gas mixing unit with flow meter and scaling in 0.5 L/min increments and overpressure valve (pressure in the tube not more than 60 mbar); flow range must be selected in accordance with the size of the Nuflow nasal cannula
- (3) Active respiratory gas humidifier
- (4) Appropriate tubing system, incl. inspiration tube (heated) with patient-side connection, ID 10 mm, and humidifier chamber


Figure 1 - Nuflow nasal cannula in combination with the medinBLENDER (REF 1085_15, REF 1090, 1091, 1093, 1094)

- **Combination with a medin nCPAP driver with High Flow mode**

- (1) Nuflow nasal cannula incl. fixation patch
- (2) nCPAP driver with High Flow mode and overpressure protection
- (3) Active respiratory gas humidifier
- (4) Appropriate tubing system, incl. inspiration tube (heated) with patient-side connection, ID 10 mm, and humidifier chamber


Figure 2 - Nuflow nasal cannula in combination with medin nCPAP driver with High Flow mode (REF 1080, 3080, 3000)

• **Combination with a ventilator**

- (1) Nuflow nasal cannula incl. fixation patch
- (2) Ventilator with High Flow mode and overpressure protection
- (3) Active respiratory gas humidifier
- (4) Appropriate tubing system, single or double tubing system (heated) with patient-side connection, ID 10 mm or ID 15 mm (Y-connection), and humidifier chamber


Figure 3 - Nuflow nasal cannula in combination with a ventilator

4.2 Structure and set-up of the system


Figure 4 - Set-up of Nuflow nasal cannula, device side

Warning:

If the packaging of a part is damaged or if the products themselves are damaged, these parts should not be used and instead should be disposed of.

Warning:

The use of chin straps or films to actively close the mouth is not permissible.

Warning:

The use of creams can have a negative effect on the fixation of the Nuflow nasal cannula.

Warning:

The Nuflow nasal cannula may be secured by a medical fixation plaster.


Warning:

Do not use the Nuflow nasal cannula if any connections are loose.

Warning:

Please check the patency of the Nuflow nasal cannula before using it on the patient by setting a flow on your suitable high-flow therapy device.

(1) Connect the tubing set to a suitable device for high-flow therapy and to the humidifier. If available, insert the temperature sensors as far as they will go; their tip must be in the gas flow.

(2) Connect the Nuflow nasal cannula to the venti lati on tube. Use the adapter provided, if necessary.

Figure 5 - Fixation patch, small

(3) Prepare the provided small fixation patch for Nuflow sizes small and medium or the large fixation patch for Nuflow sizes large and x-large. Before the fixation patch is adhered to the skin, it should be at body temperature.

(4) The skin should be dry and clean to ensure an optimal hold.


Figure 6 - Fixation patch, large

Note:

When securing the hydrocolloid fixation, ensure that it is not adhered to the eye area, in particular not to the area of the lower eyelids.

Scissors can be used to slightly adapt the hydrocolloid fixation to the patient's anatomy.

5) Select the size of the Nuflow nasal cannula such that the prongs fill a max. of 50% of the respective nostril!

Warning:

The prongs of the Nuflow nasal cannula should not fill the nostril; at least 50% of the nostril should remain clear to ensure the flushing effect and to avoid uncontrolled pressure increases which can harm the patient.


Figure 7 - Selecting the appropriate size of the Nuflow nasal cannula

(6) Release the fixation patch from the backing film. Secure this to the patient's left and right cheeks. Take note of the information on the backing of the fixation indicating the correct side (L = left ; R = right; the sides indicated are from the patient's perspective).


Figure 8 - Releasing the fixation patch from the backing film (1) and fixation to the cheek (2)

(7) There is a thin removable film on the inner side of the fixation patch. Remove this before you position the Nuflow nasal cannula.

(8) Open the tube fixation of the Nuflow nasal cannula and pull the cannula over the patient's head. Ensure that the fixation tabs which secure the Nuflow nasal cannula are open.


Figure 9 - Removing the removable film of the fixation patch

(9) Guide the prongs of the Nuflow nasal cannula into the nose. Close the upper part of the fixation patch by gently pulling downwards.

(10) Then gently push the tube fixation in the direction of the head until stable guidance of the tubes is ensured.


Figure 10 - Securing the Nuflow nasal cannula by closing the fixation patch

(11) Using gentle pressure, press the fixation tabs of the fixation bandage with your finger in the direction of the Nuflow nasal cannula to ensure optimal contact of the adhesive surface with the Nuflow nasal cannula.

(12) The fixation tabs of the fixation patches can be opened and closed multiple times. The fixation patches can remain on the patient's skin during pauses/interruptions in the high-flow therapy.

Figure 11 - Secured Nuflow nasal cannula

Fixation patches which loosen at the edges (e.g. at the nose) do not necessarily need to be replaced. If the fixation patch continues to adhere and the Nuflow nasal cannula is still secure, the fixation patch can remain on the skin.

Note:

Please ensure that the Nuflow nasal cannula does not completely slide into the nose. This can cause the nostrils to become sealed and at the same time cause redness of or even injury to the skin due to excessive pressure on the skin.

Warning:

Ensure that the weight of the ventilation tube does not pull on the Nuflow nasal cannula or the patient's head.

Warning:

The application needs to be regularly checked for the presence of secretions as well as the removal of secretions.

Warning:

Ensure that the patient does not lie on the tube.

Warning:

If the adapter included with the Nuflow nasal cannula is not used, this should not be kept within the patient's reach (small part which can be swallowed).


Figure 12 - Tube fixation

The optimal positioning of the supply tubes of the Nuflow nasal cannula is close to the patient's skin. Humidified breathing gas is heated by the patient's body temperature and thus the formation of condensate is reduced.

5 Products which can be connected to the Nuflow

The Nuflow shall be combined only with the following devices tested by medin for compatibility:

Product name	REF	Description
medinCNOmini	3080	Device to provide the set volume flow rate
medin-NC3	3000	Device to provide the set volume flow rate
medinSINDI	1080	Device to provide the set volume flow rate
medinBLENDER	1090 1091 1093 1094	Device to provide the set volume flow rate
F&P 900MR850	MR850ARU	
F&P 900MR810	MR810ARU	
10x nCPAP circuit set (single - limb)	1207MKI	With humidifier chamber
10x nCPAP circuit set (single - limb)	1207	Without humidifier chamber

The use of products that are not compatible puts the patient at risk and may impair the proper and correct functioning of the product.

6 Accessories

There are no accessories for the Nuflow.

7 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

8 Symbols

The symbols explained below are used as labels on the product or in these instructions for use:

Symbol	Meaning	Symbol	Meaning
CE	CE mark		Temperature limit
MD	Medical device		Latex-free
	Follow instructions for use		Manufacturer
REF	Article number		Protect from sunlight
LOT	Production lot number, batch		Date of manufacture
	Do not reuse		Use before
	Keep dry		Quantity
XX REP	Authorized representative		Do not use if packaging is Damaged
RI	Revision index	Rx ^{USA} _{only}	Prescription only (US)

9 Abbreviations

Abbreviation	Meaning
%	Percent
(m)bar	(Milli-) bar
(n)CPAP	(nasal) Continuous Positive Airway Pressure
°C	Degrees Celsius
G	Gram
ID	Inner diameter
L/min	Liter per minute
mm	Millimeter
OD	Outer diameter

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© 2023 medin Medical Innovations GmbH. All rights reserved. ドイツ印刷