


INDICATIONS FOR USE:

The Oxygen Nasal Cannula Loop is a single patient use oxygen delivery system consisting of a silicone nasal cannula body.



WARNINGS:

- Intended for use with instruction from qualified medical personnel.
- Do not use any hydrocarbon based substance on any part of this device.
- Do not administer supplemental oxygen in the presence of fire, sparks and open flame due to potential fire hazard.

CAUTIONS:

- Read and understand the contents of this insert and demonstrate proficiency in the application of this device prior to use.
- Do not use the product if it shows signs of damage or fails to perform as indicated.
- Single patient use; not intended to be reprocessed or reused on multiple patients as it may induce cross contamination.
- The tube should be checked periodically and replaced when damage, blockage or discoloration occurs.
- Please use this product within 5 years from the date shown on the date lot label (Lot YY-MM-DD).

INSTRUCTIONS FOR USE

1. Connect the soft connector on the tubing to the adaptor on the Ventilator/ blender/ flow driver.
2. Attach oxygen tubing to an oxygen flow meter and adjust the oxygen flow.
CAUTION: To avoid any operation error that may result in lung injury, adjust the oxygen flow before connecting to a patient and remove the nasal cannula before turning off the ventilator.
3. Recheck the gas flow system to assure it is functioning correctly and gas flow is being delivered to the Oxygen Nasal Cannula.
CAUTION: Place a hand in front of the tip to check for the air flow and other areas for leakage.
4. Gently slip the tip into the nostrils and rest the nasal cannula on the philtrum.
5. Re-adjust the position of the cannula as needed.
CAUTION: When applying with a flow rate larger than 4 LPM, humidifier should be used to minimize patient's discomfort.

SPECIFICATIONS

Connections:

ISO 10 ID soft connector

Operation Temperature: -18°C - 50°C (-0.4°F - 122°F)

Storage Temperature: -40°C - 60°C (-40°F - 140°F)

All materials meet ISO-10993 bio-compatibility requirements.

ORDER INFORMATION

REF	Name	Size	Tube Length	Flow rate
1300	Oxygen Nasal Cannula 0	0	35cm±1cm	0 – 4 lpm
1301	Oxygen Nasal Cannula 0.5	0.5	35cm±1cm	0 – 4 lpm
1302	Oxygen Nasal Cannula 1	1	35cm±1cm	0 – 7 lpm

Distributed by :

medin Medical Innovations GmbH

Adam-Geisler-Str.1, 82140 Olching, Germany

Tel.:+49 (0) 8142-448460 support@medingmbh.com



GaleMed Corporation

87, Li-gong 2nd Road, Wu-jia, I-Lan 268 Taiwan



Obelis s.a. Bd. Général Wahis 53, 1030 Brussels, Belgium