
Instructions for use

nCPAP circuit set

Contents

EN

DE

CZ

DA

EL

ES

FI

FR

HU

IT

JA

NL

NO

PL

PT

RO

RU

SL

SV

TR

ZH

The user has to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State/Country in which the user is established.

Serious incidents are defined as malfunctions, failures or changes in characteristics or performance, or inadequacies in the labelling or instructions for use of a medical device, which directly or indirectly have led, could have led or could lead to the death or serious deterioration of the state of health of a patient, user or other person.

For nCPAP circuit set

CE mark according to the current conformity assessment procedure and declaration of conformity (see website).

Classification (MDR (EU) 2017/745 Annex VIII): IIA

medin Medical Innovations GmbH

Adam-Geisler-Str. 1 - 82140 Olching - Germany

+49 8142 448 460 - www.medin-medical.com - info@medin-medical.com

Made in Germany

IFU_nCPAP_circuit_set_Rev11_EN

REF 206746, 1207, 1207MKI

© 2023 medin Medical Innovations GmbH. All rights reserved. Printed in Germany.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of medin Medical Innovations GmbH.

This document may be revised, replaced, or made obsolete by other documents by medin Medical Innovations GmbH at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the customer support department of medin Medical Innovations GmbH or download the current applicable version from our homepage www.medin-medical.com. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way the right of medin Medical Innovations GmbH to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, medin Medical Innovations GmbH has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. medin Medical Innovations GmbH sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in the user manual.

medin Medical Innovations GmbH shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-medin medical Innovations GmbH parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

medin Medical Innovations GmbH will make available, on request, circuit diagrams, component parts lists, descriptions, calibration instructions, or other information that will assist appropriately trained personnel to repair those parts of the equipment designated by medin Medical Innovations GmbH to be repairable.

For all proprietary as well as third-party trademarks used by medin Medical Innovations GmbH, see www.medin-medical.com/de/Patents.html. Product and/or company names marked with a ® may be the trademarks and/or registered trademarks of their respective owners.

Table of contents

1	Introduction.....	
1.1	About these instructions for use	
1.2	General liability conditions	
2	Intended use.....	
2.1	nCPAP circuit set.....	
2.2	Indications, expected clinical benefit, contraindications, side effects	
3	Technical specifications	
3.1	nCPAP circuit sets	
3.1.1	REF 1207 – nCPAP circuit set (single-limb)	
3.1.2	REF 1207MKI – nCPAP circuit set (single-limb).....	
3.1.3	REF 206746– nCPAP circuit set (dual-limb)	
3.2	Environmental conditions (operation / transport / storage).....	
4	Principles of function, assembly and use	
5	Products which can be connected to the nCPAP circuit sets	
6	Accessories	
7	Disposal	
8	Symbols	
9	Abbreviations	

1 Introduction

1.1 About these instructions for use

These instructions for use contain information on starting up and using the nCPAP circuit sets. They also contain safety information, describe the product function, and give an overview of the accessories needed.

Users must be very familiar with the information and warnings given in these instructions for use in order to use the nCPAP circuit sets safely. Nonetheless, they are no substitute for training. The instructions and warnings are categorized as follows:

Warning: Warnings must be followed in order to prevent possible serious consequences for the patient or the user.

Attention: This indicates hazards which could damage the device or impair its functionality.

Note: This indicates how the nCPAP circuit sets can be used more efficiently.

Keep these instructions for use close to the nCPAP circuit sets and easily accessible.

Please contact the dealer or go directly to the manufacturer if you have any questions or comments about these instructions for use.

Warning:

- Read these instructions for use carefully before using the products, and follow the instructions, warnings and tips closely.
- The nCPAP circuit sets may only be used by individuals with a thorough knowledge of and in accordance with these instructions for use for the purpose described under the section headed "Intended use".
- The use of devices and consumables used in combination with the nCPAP circuit sets is described in the respective instructions for use. Users must also comply with these instructions for use.

1.2 General liability conditions

The general terms of trade of medin Medical Innovations GmbH (hereinafter referred to as medin) are binding.

medin accepts no responsibility or liability for the safe use of the nCPAP circuit sets if an operation that does not comply with intended use is performed. The user is responsible for respecting the specifications in the instructions for use. In addition, the use of the nCPAP circuit sets requires training of the user by a qualified person.

The duration of use of the nCPAP circuit sets is defined as a maximum of 7 days.

medin accepts no liability for damage caused by failure to comply with these instructions.

2 Intended use

2.1 nCPAP circuit set

The non-invasive circuit set is a single-use product for non-invasive respiratory support and oxygen therapy from newborn infants up to young children. The circuit set transports, heated humidified gas from a humidifier to a patient interface.

Warning:

- DO NOT – reprocess this single use circuit
- DO NOT – soak, rinse, wash or sterilize this product
- DO NOT – pull, stretch or “milk” the circuit tubing
- DO NOT – use this circuit if gas temperatures in the outlet of the humidifier can exceed 154°F (68°C)
- DO NOT – use this circuit if any connections are loose
- DO NOT – use the product if the single unit packaging is damaged
- The circuits must NOT be allowed to sag between the patient and the humidifier, thereby preventing a siphon from developing
- The humidifier should be located lower than the patient
- The nCPAP circuits sets shall only be combined with the products listed
- The duration of use of the nCPAP circuit sets is defined as a maximum of 7 days
- The product must be disposed of in accordance with the manufacturer's instructions
- Only use the nCPAP circuit sets while monitoring SpO₂

2.2 Indications, expected clinical benefit, contraindications, side effects

Indications

Neonatal respiratory insufficiency.

Contraindications

Non-invasive respiratory therapy is contraindicated if certain clinical situations, anatomic malformations, or the degrees of severity of a disease necessitate intubation and mechanical ventilation or make non-invasive treatment impossible.

Benefit

Non-invasive respiratory therapy improves or maintains the gas exchange in the lungs of the patients.

Intended patient / target group

Newborn infants up to young children

Intended user

The non-invasive circuit set is intended to be used by trained healthcare professionals under the direct supervision of a licensed physician.

Intended environment


In a clinical environment while monitoring the patient's oxygen saturation.

3 Technical specifications

3.1 nCPAP circuit sets

3.1.1 REF 1207 – nCPAP circuit set (single-limb)

Technical Data

Nr.	Component	Specifications
1	Heated wire Kit	Electrical resistance: 22,5 Ω
2	Tubing	Length: 1100 mm (heated tubing part)
3	Temp. Cap	Sensor port: Chamber 7,6 mm ID, Patient 7,6 mm ID
4	L connector	-
5	Straight Connector	10M/10M
6	EVA Tubing	Length: 250 mm (unheated tubing part)
7	Straight Connector	10M/10F
8	Straight Connector	22F, 13B
9	EV Tubing	Length: 600 mm
10	Pressure Connector	-
11	Pressure line	4.1/2.7mm


Technical Data		
		Length: 1800mm
12	TPR Connector	-
13	Straight Connector	22M, 22F
14	Straight Connector	22M/15F,22F
15	Straight Connector	15M/15F,22F
16	Silicone tube	2mm ID/4mmOD Length: 45 mm

3.1.2 REF 1207MKI – nCPAP circuit set (single-limb)

Technical Data		
Nr.	Component	Specifications
1	Heated wire Kit	Electrical resistance: 22,5 Ω
2	Tubing	Length: 1100 mm (heated tubing part)
3	Temp. Cap	Sensor port: Chamber 7,6 mm ID, Patient 7,6 mm ID
4	L connector	-

Technical Data		
5	Straight Connector	10M/10M
6	EVA Tubing	Length: 250 mm (unheated tubing part)
7	Straight Connector	10M/10F
8	Straight Connector	22F, 13B
9	EV Tubing	Length: 600 mm
10	Pressure Connector	-
11	Pressure line	4.1/2.7mm, length: 1800mm
12	TPR Connector	-
13	Straight Connector	22M, 22F
14	Straight Connector	22M/15F,22F
15	Straight Connector	15M/15F,22F
16	Silicone tube	2mm ID/4mmOD, length 45 mm
17	Neonate chamber	22mmOD Volume: 53-130 ml, Operation temperature: 34 °C – 37 °C Chamber limb: 22mm ID, length: 600mm

3.1.3 REF 206746– nCPAP circuit set (dual-limb)

Technical Data		
		
Nr	Component	Specifications

Technical Data			
A	Inhalation Heated Wire tube with Y Connector	Y connector	10 OD
		Connector Cap	Luer
		T Connector	7.9F/13M
		Straight Connector	10M/13M,13B
		Straight Connector	13F,13B
		Heating wire	Electrical resistance: 22,5 Ω
		L connector	22ID
		Tubing	1100 mm
		Tubing	300 mm
		EVA smoothbore	300 mm
B	Exhalation non-Heated tube	Tubing	1000 mm
		Tubing	500 mm
		Straight Connector	10M/13M,13B
		Straight Connector	10F/10F
		Straight Connector	15M,13B
C	Chamber	-	
D	Water trap	Y connector 10OD	
E	Chamber limb	Straight Connector	150D,150D
		Straight Connector	150D,22ID
		Corrugate Tube	22ID, 22ID connector
F	Connector Kit	a. Straight connector	22OD, 22ID
		b. Silicone tube	3.0mmID/5.0mmOD 45mm
		c. T connector	10ID,4.3ID,10OD/7.6ID
		d. TPR connector	9.8ID,9.8ID, 2 pcs
		e. Straight connector	22/15, 15OD 2pcs

Technical Data			
		f. Straight connector	22/15, 22OD
G	Pressure line		length: 1800mm optional

3.2 Environmental conditions (operation / transport / storage)

Environmental conditions for the nCPAP circuit sets during operation:


Temperature:	15°C to 40°C
Cleanliness:	The hygiene regulations valid for the hospital are to be observed. In addition, the nCPAP circuit sets are to be used under standard hospital ambient conditions.

Environmental conditions for the nCPAP circuit sets during storage/transport:

Temperature:	-30°C to 70°C
Cleanliness:	The nCPAP circuit sets must be protected from contamination during transport and storage.
Other:	The nCPAP circuit sets generator must be stored dry and protected from sunlight.

4 Principles of function, assembly and use

EN

1. Insert the chamber into the Humidifier.
2. Connect tube b to the ventilator and the chamber. If necessary use the additional connectors d.
3. Connect tube a to the chamber.
4. Plug in the cable for the heating wire.
5. Insert the two temperature sensors. Make sure they are inserted far enough – their tips have to reach into the gas flow.
6. Attach the Medijet/ Miniflow generator to tube a.
7. MINIFLOW: If a pressure measurement is necessary: add the t-piece e between Miniflow and the inspiration tube and connect the pressure line c with the T-piece and the ventilator.
8. MEDIJET: Connect the pressure line c with the driver and the Medijet.
9. Miniflow: Connect Miniflow and the ventilator with tube d.
10. Test the function and the tightness without contact to the patient according to Medijet's/ Miniflow's user manual.
11. The system is ready for use.
12. The system is ready for use and safe if no alarms are present.

5 Products which can be connected to the nCPAP circuit sets

The nCPAP circuit sets shall be combined only with the following devices tested by medin for compatibility:

Table 1 - Combinable products

Product name	REF	Combinations
20x Medijet 10x Medijet	1000 10162670	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
20x Miniflow 10x Miniflow	4000 10162671	206746 - 10x nCPAP circuit set (dual-limb)
10x Nuflow nasal cannula (small) 10x Nuflow nasal cannula (medium) 10x Nuflow nasal cannula (large) 10x Nuflow nasal cannula (x-large) Nuflow starter set s/m Nuflow starter set l/xl	1320 1330 1340 1350 1393 1394	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinCNO medinCNOMini	3090 3080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medin-NC3	3000	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinSINDI	1080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinBLENDER	1090 1091 1093 1094	206746 - 10x nCPAP circuit set (dual-limb)
Bubble Water Column with overpressure valve	2040	206746 - 10x nCPAP circuit set (dual-limb)
F&P 900MR850	MR850ARU	All nCPAP circuit sets
F&P 900MR810	MR810ARU	All nCPAP circuit sets

6 Accessories

There are no accessories for the nCPAP circuit sets.

7 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

8 Symbols

The symbols explained below are used as labels on the product or in these instructions for use:

Table 2 - Symbols

Symbol	Meaning	Symbol	Meaning
	CE mark		Temperature limit
	Medical device		Latex-free
	Follow instructions for use		Manufacturer
	Article number		Protect from sunlight
	Production lot number, batch		Date of manufacture
	Do not reuse		Use before
	Keep dry		Quantity
	Authorized representative		Do not use if packaging is Damaged
	Revision index		

9 Abbreviations

Table 3 - Abbreviations

Abbreviation	Meaning
%	Percent
(n)CPAP	(nasal) Continuous Positive Airway Pressure
°C	Degrees Celsius
ID	Inner diameter
mm	Millimeter
OD	Outer diameter
Ω	Ohm

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Tabelle 1 - Kombinierbare Produkte

Produktbezeichnung	REF	Kombinationen
20x Medijet 10x Medijet	1000 10162670	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
20x Miniflow 10x Miniflow	4000 10162671	206746 - 10x nCPAP circuit set (dual-limb)
10x nasal cannula (small) 10x nasal cannula (medium) 10x nasal cannula (large) 10x nasal cannula (x-large) Nuflow starter set s/m Nuflow starter set l/xl	1320 1330 1340 1350 1393 1394	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinCNO medinCNOMini	3090 3080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medin-NC3	3000	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinSINDI	1080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinBLENDER	1090 1091 1093 1094	206746 - 10x nCPAP circuit set (dual-limb)
Bubble water column with overpressure valve	2040	206746 - 10x nCPAP circuit set (dual-limb)
F&P 900MR850	MR850ARU	Alle nCPAP circuit sets
F&P 900MR810	MR810ARU	Alle nCPAP circuit sets

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Tabella 1 - Prodotti combinabili

Nome del prodotto	REF	Combinazioni
20x Medijet 10x Medijet	1000 10162670	207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
20x Miniflow 10x Miniflow	4000 10162671	206746 - 10x nCPAP circuit set (dual-limb)
10x Nuflow nasal cannula (small) 10x Nuflow nasal cannula (medium) 10x Nuflow nasal cannula (large) 10x Nuflow nasal cannula (x-large) Nuflow starter set s/m Nuflow starter set l/xl	1320 1330 1340 1350 1393 1394	207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinCNO medinCNOMini	3090 3080	207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medin-NC3	3000	207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinSINDI	1080	207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinBLENDER	1090 1091 1093 1094	206746 - 10x nCPAP circuit set (dual-limb)
Bubble water column with overpressure valve	2040	206746 - 10x nCPAP circuit set (dual-limb)
F&P 900MR850	MR850ARU	Tutti i set circuito nCPAP
F&P 900MR810	MR810ARU	Tutti i set circuito nCPAP

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製品名	REF	組み合わせ
20x Medijet 10x Medijet	1000 10162670	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
20x Miniflow 10x Miniflow	4000 10162671	206746 - 10x nCPAP circuit set (dual-limb)
10x Nuflow nasal cannula (small) 10x Nuflow nasal cannula (medium) 10x Nuflow nasal cannula (large) 10x Nuflow nasal cannula (x-large) Nuflow starter set s/m Nuflow starter set l/xl	1320 1330 1340 1350 1393 1394	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinCNO medinCNOmini	3090 3080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medin-NC3	3000	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinSINDI	1080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinBLENDER	1090 1091 1093 1094	206746 - 10x nCPAP circuit set (dual-limb)
Bubble water column with overpressure valve	2040	206746 - 10x nCPAP circuit set (dual-limb)
F&P 900MR850	MR850ARU	すべての nCPAP 回路セット
F&P 900MR810	MR810ARU	すべての nCPAP 回路セット

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
3.1.3 REF 206746– nCPAP circuit set (dual-limb)

Technische gegevens		
 The diagram illustrates the nCPAP circuit set (dual-limb). It includes a main circuit assembly (A) with a long tube (B) and a humidifier (C). The air inlet (D) is connected to the main assembly. The circuit is terminated with connectors (E and F). A long tube (G) is also shown. Below the diagram, small images of components are labeled a through f.		
Nr.	Onderdeel	Specificaties

Technische gegevens			
A	Verwarmde inhalatiedraadslang met Y-stuk connector	Y-stuk connector	10 OD
		Connectordop	Luer
		T-stuk connector	7.9F/13M
		Rechte connector	10M/13M,13B
		Rechte connector	13F,13B
		Verwarmingsdraad	Elektrische weerstand: 22,5 Ω
		L-connector	22ID
		Slang	1100 mm
		Slang	300 mm
		EVA smoothbore slang	300 mm
B	Niet-verwarmde uitademingsslang	Slang	1000 mm
		Slang	500 mm
		Rechte connector	10M/13M,13B
		Rechte connector	10F/10F
		Rechte connector	15M,13B
C	Kamer		-
D	Waterslot		Y-stuk connector 10OD
E	Kamerslang	Rechte connector	15OD,15OD
		Rechte connector	15OD, 22ID
		Ribbelslang	22ID, 22ID connector
F	Aansluitkit	a. Rechte connector	22OD, 22ID
		b. Siliconen slang	3,0 mm ID/5,0 mm OD 45 mm
		c. T-stuk connector	10ID,4.3ID,10OD/7.6ID
		d. TPR-connector	9.8ID,9.8ID, 2 stuks
		e. Rechte connector	22/15, 15OD 2 stuks

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
Tabel 1 - combineerbare producten

Productnaam	REF	Combinaties
20x Medijet 10x Medijet	1000 10162670	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
20x Miniflow 10x Miniflow	4000 10162671	206746 - 10x nCPAP circuit set (dual-limb)
10x Nuflow nasal cannula (small) 10x Nuflow nasal cannula (medium) 10x Nuflow nasal cannula (large) 10x Nuflow nasal cannula (x-large) Nuflow starter set s/m Nuflow starter set l/xl	1320 1330 1340 1350 1393 1394	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinCNO medinCNOmini	3090 3080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medin-NC3	3000	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinSINDI	1080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinBLENDER	1090 1091 1093 1094	206746 - 10x nCPAP circuit set (dual-limb)
Bubble water column with overpressure valve	2040	206746 - 10x nCPAP circuit set (dual-limb)
F&P 900MR850	MR850ARU	Alle nCPAP-circuitsets
F&P 900MR810	MR810ARU	Alle nCPAP-circuitsets

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Tabel 1 - produse care pot fi combinate

Denumirea produsului	REF	Combinații
20x Medijet 10x Medijet	1000 10162670	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
20x Miniflow 10x Miniflow	4000 10162671	206746 - 10x nCPAP circuit set (dual-limb)
10x Nuflow nasal cannula (small) 10x Nuflow nasal cannula (medium) 10x Nuflow nasal cannula (large) 10x Nuflow nasal cannula (x-large) Nuflow starter set s/m Nuflow starter set l/xl	1320 1330 1340 1350 1393 1394	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinCNO medinCNOMini	3090 3080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medin-NC3	3000	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinSINDI	1080	1207 - 10x nCPAP circuit set (single-limb) 1207MKI - 10x nCPAP circuit set (single-limb)
medinBLENDER	1090 1091 1093 1094	206746 - 10x nCPAP circuit set (dual-limb)
Bubble water column with overpressure valve	2040	206746 - 10x nCPAP circuit set (dual-limb)
F&P 900MR850	MR850ARU	Toate seturile de circuite nCPAP
F&P 900MR810	MR810ARU	Toate seturile de circuite nCPAP

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© 2023 medin Medical Innovations GmbH. All rights reserved. ドイツ印刷