
Instructions for use

Medijet

Contents

EN

DE

CZ

DA

EL

ES

FI

FR

HU

IT

JA

NL

NO

PL

PT

RO

RU

SL

SV

TR

The user has to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State/Country in which the user is established.

Serious incidents are defined as malfunctions, failures or changes in characteristics or performance, or inadequacies in the labelling or instructions for use of a medical device, which directly or indirectly have led, could have led or could lead to the death or serious deterioration of the state of health of a patient, user or other person.

For Masks:

CE mark according to the current conformity assessment procedure and declaration of conformity (see website).

Classification (MDR (EU) 2017/745 Annex VIII): I

For Medijet and Prongs:

CE mark according to the current conformity assessment procedure and declaration of conformity (see website).

Classification (MDR (EU) 2017/745 Annex VIII): IIA

medin Medical Innovations GmbH
Adam-Geisler-Str. 1 - 82140 Olching - Germany
+49 8142 448 460 - www.medin-medical.com - info@medin-medical.com
Made in Germany

IFU_Medijet_Rev17_EN

Medijet:

REF 1000, 10162670, 1000-10

Prongs:

REF 1200-01, 1200-02, 1200-03, 1200-21, 1200-22, 1200-32, 1200-33, 10162662

Masks:

REF 1200-04, 1200-05, 1200-06, 1200-07, 1200-08, 1200-09, 10162666

© 2023 medin Medical Innovations GmbH. All rights reserved. Printed in Germany.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of medin Medical Innovations GmbH.

This document may be revised, replaced, or made obsolete by other documents by medin Medical Innovations GmbH at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the customer support department of medin Medical Innovations GmbH or download the current applicable version from our homepage www.medin-medical.com. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way the right of medin Medical Innovations GmbH to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, medin Medical Innovations GmbH has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. medin Medical Innovations GmbH sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in the user manual.

medin Medical Innovations GmbH shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-medin medical Innovations GmbH parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

medin Medical Innovations GmbH will make available, on request, circuit diagrams, component parts lists, descriptions, calibration instructions, or other information that will assist appropriately trained personnel to repair those parts of the equipment designated by medin Medical Innovations GmbH to be repairable.

For all proprietary as well as third-party trademarks used by medin Medical Innovations GmbH, see www.medin-medical.com/de/Patents.html. Product and/or company names marked with a [®] may be the trademarks and/or registered trademarks of their respective owners.

Table of contents

1	Introduction.....	
1.1	About these instructions for use	
1.2	General liability conditions	
2	Intended use.....	
2.1	Medijet	
2.2	Prong	
2.3	Mask	
2.4	Indications, expected clinical benefit, contraindications	
3	Principles of function and technical specifications.....	
3.1	Medijet	
3.2	Bonnets, prongs, masks, Medifoam	
3.2.1	Bonnets.....	
3.2.2	Prongs.....	
3.2.3	Masks.....	
3.2.4	Medifoam	
3.3	Environmental conditions (operation / transport / storage).....	
4	Assembly and use	
4.1	Equipment needed	
4.2	Additional aids.....	
4.3	Structure and set-up of the system	
5	Products which can be connected to the Medijet.....	
5.1	Medijet	
5.2	Prongs, masks	
6	Accessories	
7	Disposal	
8	Symbols	
9	Abbreviations	

1 Introduction

1.1 About these instructions for use

These instructions for use contain information on starting up and using the Medijet/prong/mask. They also contain safety information, describe the product function, and give an overview of the accessories needed.

Users must be very familiar with the information and warnings given in these instructions for use in order to use the Medijet/prong/mask safely. Nonetheless, they are no substitute for training. The instructions and warnings are categorized as follows:

Warning: Warnings must be followed in order to prevent possible serious consequences for the patient or the user.

Attention: This indicates hazards which could damage the device or impair its functionality.

Note: This indicates how the Medijet/prong/mask can be used more efficiently.

Keep these instructions for use close to the Medijet/prong/mask and easily accessible.

Please contact the dealer or go directly to the manufacturer if you have any questions or comments about these instructions for use.

Warning:

- Read these instructions for use carefully before using the products, and follow the instructions, warnings and tips closely.
- The Medijet/prong/mask may only be used by individuals with a thorough knowledge of and in accordance with these instructions for use for the purpose described under the section headed "Intended use".
- The use of devices and consumables used in combination with the Medijet/prong/mask is described in the respective instructions for use. Users must also comply with these instructions for use.

1.2 General liability conditions

The general terms of trade of medin Medical Innovations GmbH (hereinafter referred to as medin) are binding.

medin accepts no responsibility or liability for the safe use of the Medijet/prong/mask if an operation that does not comply with intended use is performed. The user is responsible for respecting the specifications in the instructions for use. In addition, the use of the Medijet/prong/mask requires training of the user by a qualified person.

The duration of use of the Medijet/prong/mask is defined as a maximum of 7 days.

medin accepts no liability for damage caused by failure to comply with these instructions.

2 Intended use

2.1 Medijet

The medijet is a single-use product for non - invasive respiratory support and oxygen therapy of newborn infants. It is used in combination with a CPAP driver which supplies the necessary gas flow to generate the CPAP pressure. The medijet is used under the supervision of competent personnel trained for this purpose in a clinical environment.

Warning:

- The oxygen saturation in the patient's blood must be continuously monitored during use.
- The Medijet, masks, prongs, bonnets and Medifoam are disposable products; they may be used only on a single patient – and may not be reprocessed.
- The Medijet may only be used in combination with prongs, masks, bonnets and Medifoam
- The duration of use of the Medijet/prong/mask is defined as a maximum of 7 days.
- The product must be disposed of in accordance with the manufacturer's instructions

2.2 Prong

The prong is a single-use product for non-invasive respiratory support and oxygen therapy of newborn infants. It is used in combination with a medijet or miniflow. The prong is used under the supervision of competent personnel trained for this purpose in a clinical environment.

2.3 Mask

The mask is a single-use product for non - invasive respiratory support and oxygen therapy of newborn infants. It is used in combination with a medijet or miniflow. The mask is used under the supervision of competent personnel trained for this purpose in a clinical environment.

2.4 Indications, expected clinical benefit, contraindications

Indications

Neonatal respiratory insufficiency

Contraindications

Non-invasive respiratory support is contraindicated if certain clinical situations, anatomic malformations, or the degrees of severity of a disease necessitates intubation and mechanical ventilation or make non-invasive treatment impossible

Benefit

The medijet/prong/mask improves or maintains the gas exchange in the lungs of the patients

Intended patient / target group

Neonatal patients

Intended user

Medijet/prong/mask is intended to be use by trained health care professionals under the direct supervision of a licensed physician

Intended environment

In a clinical environment

3 Principles of function and technical specifications

3.1 Medijet

- The Medijet is a patient interface which can be used in combination with an nCPAP driver approved by medin.
- The Medijet can be connected via a single tubing system to an nCPAP driver or ventilator. The CPAP parameters are adjusted on the nCPAP driver or ventilator. The nCPAP driver or the ventilator provides the necessary gas flow. The pressure meter tube must be connected to the pressure meter connection of the nCPAP driver or the ventilator. The Medijet offers a medication nebulizer connection as an option.
- The prong adapter of the Medijet is tiltable. This allows the direction of the prong adapter to be adjusted to the baby's nasal anatomy and weight which facilitates fixation. A prong or a mask selected according to the baby's size is placed on the prong adapter. The prong or the mask ensures a soft contact between the Medijet and the baby's nose and can be used in alternation to avoid pressure points.
- The Medijet itself is secured to the bonnet which the baby wears.

Technical data	
General	
Connection size	M10/ F7.4
Weight	10 g
Pressure connection	F4,3
Duration of use	Max. 7 days
Parameters & sensors	
CPAP pressure	0 to > 20 mbar
Principles of function	Benveniste valve
Noise level	61 dBA
Volume reservoir	4 ml - constant flow
nCPAP pressure	≥ 4.5 mbar at 7 l/min flow
Connections	
Inspiration	OD: 9.9 – 10.9 mm, conical / ID: 7.5 mm
Pressure line	Luer-like port

3.2 Bonnets, prongs, masks, Medifoam

3.2.1 Bonnets

A divided fleece is sewn onto the upper side of the bonnets. A fixation pillow (Medifoam) is placed between the fleece parts in which the Medijet is placed and secured with a separate Velcro strap. The straps additionally included are used to stabilize the mask or prong position and, with their Velcro ends, they can be placed anywhere on the bonnet fabric. The masks and prongs each contain loops for the fixation straps of the bonnet. In addition, they contain a pair of holes for the connection to the Medijet in which the prong adapter of the Medijet is inserted.

3.2.2 Prongs

The binasal prongs of medin (see Figure 1) serve as an interface between the Medijet and the patient's nasal area and thus complete the nCPAP system. They are introduced into the nostrils and thus establish the connection to the system. The prongs are available in seven sizes: XS, S, M, MW, L, LW, XL. Two short tubes are attached to a base piece. On the side of the base piece are fixation tabs through which the fixation straps of the bonnets are pulled in order to secure the prong to the bonnet. The base of the tubes has a tapered shape in order to seal the transition between the nostril and the prong through more flexible adjustment.


Figure 1 - Prong with fixation tabs (1) and tapered tubes (2)

The various prong sizes differ in diameter and distance between the tubes. The appropriate prong can be selected for the patient using a special measuring tape, in order to keep skin damage on the patient to a minimum. The tabs of the prong are attached at a slanted angle. This ensures a comfortable position on the patient's cheek during fixation to the bonnet and additionally prevents the tubes from spreading when the tabs are firmly pulled apart. The prongs can be used on the Medijet and also on the Miniflow.

3.2.3 Masks

Binasal masks (see Figure 2), like the binasal prongs, serve as an interface between the Medijet and the patient's nose. medin offers masks in five sizes: XS, S, M, L and XL. They are placed on the patient's nose and are intended to form a tight seal. The masks have a malleable cushion which establishes contact with the skin. This cushion can compensate for uneven areas and because of its malleability, it reduces skin damage. The fixation straps are secured in the tabs attached to the side and the straps are then secured to the bonnet. The masks can be used on the Medijet and also on the Miniflow.


Figure 2 - Mask with connection opening (1) for the nCPAP interface, fixation tabs (2) and application cushion (3)

3.2.4 Medifoam

The Medifoam is a fixation pillow into which the Medijet is fitted. The opening in the Medifoam corresponds to the diameter of the Medijet tube. The Medifoam is secured to the bonnet together with the Medijet (see Figure 3). The Velcro strap corresponds to the dimensions of the Velcro pads on the bonnet.

Figure 3 - Medifoam (1) and Medijet fitted into Medifoam (2)

3.3 Environmental conditions (operation / transport / storage)

Environmental conditions for Medijet, prong, mask during operation:

Temperature:	15°C to 40°C
Cleanliness:	The hygiene regulations valid for the hospital are to be observed. In addition, the Medijet/prong/mask is to be used under standard hospital ambient conditions

Environmental conditions for Medijet, prong, mask during storage/transport:

Temperature:	-30°C to 70°C
Cleanliness:	The Medijet/prong/mask must be protected from contamination during transport and storage
Other:	The Medijet/prong/mask must be stored dry and protected from sunlight

4 Assembly and use

4.1 Equipment needed

To perform effective CPAP therapy with the Medijet, the following equipment is needed:

- An nCPAP driver (flow-controlled) or a special ventilator which provides a connection for an external measurement of the CPAP pressure and which is suitable for operation with the Medijet.
- An active respiratory gas humidifier
- A Medijet
- A Medifoam
- An appropriate tubing set which must contain an inspiration tube and a pressure measurement line. In addition, a chamber for the humidifier is needed. The connection on the tubing end of the inspiration tube must be F10/M7.4 and the pressure measurement line must contain a Luer plug (M4.3) to fit on the Medijet.
- A bonnet – its size must be selected such that it fits the baby.
- A mask or a prong which fits the size of the patient's nose.

4.2 Additional aids

- A measuring tape to help select the sizes of bonnets, masks and prongs.

Warning:

If the packaging of a part is damaged or if the products themselves are damaged, these parts should not be used and should be disposed of.

If the openings are not patent, the products in question should not be used.

4.3 Structure and set-up of the system


Figure 4 - Set-up of the Medijet


Figure 5 - Set-up of the Medijet, device side

(1) Connect the tubing set to the nCPAP driver or to the ventilator and the humidifier. If necessary, use the enclosed adapter. Insert the temperature sensors as far as they will go; their tip must be in the gas flow.

(2) Connect the Medijet to the inspiration tube (larger connector of the Medijet). Insert the pressure measurement line in the pressure measurement connector of the Medijet (smaller connector) and connect the pressure measurement line with the pressure measurement connector of the nCPAP driver or the ventilator.

Warning:

This medical device uses a small diameter connector with an alternative design that is different from the designs specified in the ISO 80369 series of standards. There is a possibility of misconnection between this medical device and a medical

device using a small diameter connector with a different alternative design, which may result in a hazardous situation and cause harm to the patient. The user must take special precautions to mitigate these reasonably foreseeable risks.

(3) Close the prong adapter of the Medijet manually and adjust the CPAP parameters, such as the pressure, on the nCPAP driver or ventilator. If this is not possible, check the tubing connections for leaks. If all of the connections are correct and the Medijet is still causing problems, it should not be used and must be disposed of.

(4) Select a prong or a mask, using the measuring tape as an aid, if necessary (see Figure 6). The prong should be selected such that it is large enough to seal the entire nostril. It is placed on the Medijet such that the arched side of the prong points toward the patient's upper lip.

Note:

The mask size should be selected such that it fits precisely along the sides of the nostrils. The mask is attached such that it conforms to the shape of the nose.

Warning:

In addition, before use, check the holes of the Medijet and the prongs/masks to ensure they are patent

Note:

Always use the largest possible prong to avoid leaks and ensure the best possible positioning – the smaller the leak is, the higher and more stable the CPAP pressure is.

Warning:

Always insert the Medijet fully into the prongs or masks and do not use any cream near the prongs, masks or Medijet since otherwise the prong or mask can slip from the Medijet.

Attention:

The prongs and masks should not be reprocessed or reused. They are only for single use. Reprocessed prongs and masks may harden and injure the patient.

(5) Select a bonnet by means of the color coding and use the measuring tape as an aid, if necessary (see Figure 6). The bonnet should be secured firmly enough but without leading to deformation of the head. Remove the separate Velcro strap from the bonnet and put it aside, within easy reach. The bonnet is placed on the head such that the gap between the two fleece parts is centered over the forehead, without exerting too much pressure on the head. The bonnet should cover the ears and be pulled down to the nape of the neck (see Figure 7). Guide the fixation straps through the loops on the prong or on the mask.


Figure 6 - Determining the size of the head (a) or the prong/mask size (b) using the measuring tape

Figure 7 - Positioning the bonnet on the patient's head

Note:

Pull the bonnet down to just above the patient's eyebrows. Bonnets which are too small lead to deformation of the patient's head. Bonnets which are too large cannot ensure sufficient fixation of the CPAP generator.

(6) Select the most appropriate angle for the prong adapter of the Medijet (45° or 60°), see Figure 8. If a prong is used, the narrower angle is generally more suitable; if a mask is used, the wider angle. Place the Medifoam fixation pillow between the two fleece pieces in the middle of the bonnet (see Figure 9). Fit the Medijet into the round opening of the foam provided. Insert the prong into the nostrils. In doing so, maintain a small distance between the nasal septum and the base of the prong. When a mask is used, ensure that it fits evenly on the skin and that a distance is maintained between the inner rim of the mask and the nostrils. Using the Velcro strap, secure the Medijet from above to the fleece pieces of the bonnet. Secure the Velcro area of the straps in the area of the ears as an extension of the fixation tabs of the prong. Do not pull the straps too tightly; they should only provide the lateral support needed.

Figure 8 - Adjusting the angle of the Medijet

Figure 9 - Securing the Medijet to the bonnet

Attention:

Ensure that you mount mask/prong in the correct way, see figure 8

Attention:

- Do not insert the prongs fully into the nostrils. The connector block should not touch the nose.
- Regularly switching between prong and mask is an effective method for avoiding or reducing damage to the patient's skin.
- The Medijet is held primarily by means of the bonnet and the foam. If the leak is too large, check the selected prong or mask size – pulling the straps more tightly does not by itself solve the problem.
- Do not pull the straps too tight. Excessively tightening the straps or inserting the prongs too far can lead to pressure points on the baby's nose or skin.

Attention:

- Ensure that the weight of the tubes does not pull directly on the patient's head or on the interface. Use cushioning to support the tubes.

(7) Regularly check the position and fixation of the Medijet and the accessories.

Attention:

Regularly check the positioning and fixation of the Medijet and the accessories since pressure points may otherwise develop. Check the prongs and masks regularly for collected secretions since there otherwise may be an obstruction.

5 Products which can be connected to the Medijet

5.1 Medijet

The Medijet shall only be combined with the following devices tested by medin for compatibility:

Product name	REF	Description
medinCNO	3090	To provide the volume flow rate
medinCNOmini	3080	
medin-NC3	3000	To provide the volume flow rate
medinSINDI	1080	To provide the volume flow rate
10x nCPAP circuit set (single-limb)	1207 1207MKI	For connection with the driver device
F&P 900MR850	MR850ARU	Humidifier Fisher&Paykel

5.2 Prongs, masks

Prongs and masks shall only be combined with the Medijet or the Miniflow.

6 Accessories

medin provides various components as accessories for the use of the Medijet:

Product name	REF	Comments
10x Prong (x-small)	1200-01	
10x Prong (medium)	1200-02	
10x Prong (x-large)	1200-03	
10x Prong (small)	1200-21	
10x Prong (large)	1200-22	
10x Prong (medium-wide)	1200-32	
10x Prong (large-wide)	1200-33	
Medin Mixed Box Prong (1XS, 2S, 2M, 1MW, 2L, 1LW, 1XL)	10162662	
10x Mask (small)	1200-04	
10x Mask (medium)	1200-05	
10x Mask (large)	1200-06	
5x Mask (x-large)	1200-07	
10x Mask (x-small)	1200-08	
10x Mask (large +)	1200-09	
Medin Mixed Mask (1XS, 1S, 1M, 1L, 1XL)	10162666	
10x Bonnet (xx-small)	1213-10	light green
10x Bonnet (x-small)	1214-10	white
10x Bonnet (small)	1215-10	yellow
10x Bonnet (medium)	1216-10	red
10x Bonnet (large)	1217-10	light blue
10x Bonnet (x-large)	1218-10	orange
10x Bonnet (xx-large)	1219-10	light green
10x Bonnet (xxx-large)	1220-10	white
10x Medifoam	1030	Fixation cushion for Medijet

Warning:

The Medijet only functions in combination with these masks, prongs, bonnets and the Medifoam, otherwise proper function cannot be guaranteed.

7 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

8 Symbols

The symbols explained below are used as labels on the product or in these instructions for use:

Symbol	Meaning	Symbol	Meaning
	CE mark		Temperature limit
	Medical device		Latex-free
	Follow instructions for use		Manufacturer
	Article number		Protect from sunlight
	Production lot number, batch		Date of manufacture
	Do not reuse		Use before
	Keep dry		Quantity
	For use by medical prescription only		Do not use if packaging is Damaged
	Authorized representative		Revision index
	Distributed by		

9 Abbreviations

Abbreviation	Meaning
%	Percent
(m)bar	(Milli-)bar
(n)CPAP	(nasal) Continuous Positive Airway Pressure
°C	Degrees Celsius
dBA	Decibel (unit of measurement for sound volume)
g	Gram
ID	Inner diameter
L/min	Liters per minute (unit of measurement for flow)
mm	Millimeter
OD	Outer diameter

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© 2023 medin Medical Innovations GmbH. All rights reserved. ドイツ印刷