


# REF 206746, REF 206748

Circuit set for ventilator with inspiration heating wire, water trap and chamber (only REF 206746)

clean - non steril



## Variations

	Inspiration	Expiration	Chamber for humidifier
REF 206746	Heated	Non heated and incl. water trap	yes
REF 206748	Heated	Non heated and incl. water trap	no

## Specification

- Total length:
  - Humidifier to Miniflow® (A): 1550mm +/- 100mm
  - Expiration (B): 1550mm +/- 100mm
  - Pressure line (C): 1800mm +/- 50mm
  - Gas inlet to humidifier (D): 600mm +/- 50mm
  - Tubes 13mm ID
- Circuit humidifier to Miniflow (A):
  - Connector Miniflow side: 10mm ID - female
  - Connector humidifier side: 22mm ID - female
  - Material tube: PP
  - Material connectors: PP, K-Resin
- Expiration tube (B)
  - Connector Miniflow side: 10mm ID - female
  - Connector ventilator side: 22mm ID - female
  - Material tube: EVA & LDPE
  - Material connectors: PP, K-Resin
  - Water trap: PP
- Pressure line
  - Connectors both sides: 4mm OD - male
  - Material tube: PP
  - Material connectors: K-Resin
- Circuit gas supply inlet to humidifier (D)
  - Connectors both sides: 22mm ID - female
  - Material tube: PP
  - Material connectors: K-Resin
- Additional T-Piece (E):
  - 10mm ID - female; 10mm OD - male and 4,3mm ID - female
  - Material: K-Resin
- Additional connectors (F):
  - 22mm OD - male and 15mm OD - male
  - Material: PP
- Tubus connector (G):
  - Material: PP

- Heating:
  - The circuit humidifier to Miniflow (A) includes a heating wire.
  - Electric resistance 22,5 ohm
  - Electric connector: Cloverleaf shape
  - Material: PP
  - Two ports for the temperature sensors are included in the circuit (Material connectors: PS, PP)
- Adapters:
  - 22/22 (OD/ID; OD/OD)
  - Material: PP
  - Soft connectors
  - Material: TPR
  - Si- Adapter
  - Material: Silicone
- Humidification chamber (only REF 206746)
  - Volume: 53 ~ 130ml
  - Connectors: Inlet/ Outlet 22mm OD - female
  - Operation Temperature: 34°C - 37°C
  - Length tubing: 970mm
  - Material:
    - Bottle: PC
    - Tray: Aluminium
    - Tubing: PP
    - Needle: ABS
    - Cap: PE
    - Buoy: Silicone and PP
    - Ring: TPR



medin Medical Innovations GmbH

Adam-Geisler-Straße 1

82140 Olching, Germany

+49 (0)8142 448 46 -0

info@medin-medical.com

[www.medin-medical.com](http://www.medin-medical.com)



REV07  
14.08.2019

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. For all proprietary trademarks and third-party trademarks used by medin Medical Innovations GmbH see [www.medin-medical.com/trademarks](http://www.medin-medical.com/trademarks).  
© 2019 medin Medical Innovations GmbH. All rights reserved.

20674X

