



206746 Circuit set

Technical specifications

For Miniflow with ventilator or BubbleCPAP.
With inspiration heating wire, water trap and chamber.

Heated inspiration - unheated expiration

clean - non steril



Specifications

Total length:

- Humidifier to Miniflow (A): 1550mm +/- 100mm
- Expiration (B): 1550mm +/- 100mm
- Pressure line (C): 1800mm +/- 50mm
- Gas inlet to humidifier (D): 600mm +/- 50mm
- Tubes 13mm ID

Circuit humidifier to Miniflow (A):

- Connector Miniflow side: 10mm ID - female
- Connector humidifier side: 22mm ID - female
- Material tube: PP
- Material connectors: PP, K-Resin

Expiration tube (B)

- Connector Miniflow side: 10mm ID - female
- Connector ventilator side: 22mm ID - female
- Material tube: EVA & LDPE
- Material connectors: PP, K-Resin
- Water trap: PP

Pressure line

- Connectors both sides: 4mm OD - male
- Material tube: PP
- Material connectors: K-Resin

Circuit gas supply inlet to humidifier (D)

- Connectors both sides: 22mm ID - female
- Material tube: PP
- Material connectors: K-Resin

Additional T-Piece (E):

- 10mm ID - female; 10mm OD - male and 4,3mm ID - female
- Material: K-Resin

Additional connectors (F):

- 22mm OD - male and 15mm OD - male
- Material: PP

Tubus connector (G):

- Material: PP

Heating:

- The circuit humidifier to Miniflow (A) includes a heating wire.
- Electric resistance 22,5 ohm
- Electric connector: Cloverleaf shape
- Material: PP
- Two ports for the temperature sensors are included in the circuit (Material connectors: PS, PP)

Adapters:

- 22/22 (OD/ID; OD/OD)
- Material: PP
- Soft connectors
- Material: TPR
- Si- Adapter
- Material: Silicone

Humidification chamber (only REF 206746)

- Volume: 53 ~ 130ml
- Connectors: Inlet/ Outlet 22mm OD - female
- Operation Temperature: 34°C - 37°C
- Length tubing: 970mm
- Material:
 - Bottle: PC
 - Tray: Aluminium
 - Tubing: PP
 - Needle: ABS
 - Cap: PE
 - Buoy: Silicone and PP
 - Ring: TPR



medin Medical Innovations GmbH

Adam-Geisler-Straße 1

82140 Olching, Germany

☎ +49 (0)8142 448 46 -0

info@medin-medical.com

www.medin-medical.com

REV08
19.05.2025

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. For all proprietary trademarks (®) and third-party trademarks (®) used by medin Medical Innovations GmbH see www.medin-medical.com/trademarks. © 2025 medin Medical Innovations GmbH. All rights reserved.