


Miniflow

Passive/variable flow generator



- A. Inspiration
- B. Expiration
- C. Separating wall
- D. Adjustable prong adapter

Parameters and sensors

CPAP pressures	0 to > 20mbar
Noise level	41dB(A)
Volume reservoir ml	4 ml -constant flush

Connectors

Inspiration	10mm OD/7.5 ID
Expiration	10mm OD/7.5 ID

Material


Terlux, TPE-S, Makrolon

Consumables

Interfaces	Prongs and masks in various sizes
Bonnets	Color-coded and different sizes available

Combinable devices

medinBLENDER



Material Miniflow:
Terlux (1), TPE-S (2), Makrolon (3)

Material Measuring Tape:
ECF Paper

Material Minifoam:
Polyurethan-Etherfoam

Material Blister:
PE, Paper

Material Packaging:
Carton



medin Medical Innovations GmbH
Adam-Geisler-Straße 1
82140 Olching, Germany
☎ +49 (0)8142 448 46 -0
✉ info@medin-medical.com
www.medin-medical.com

REV02
27.05.2025

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. For all proprietary trademarks (®) and third-party trademarks (®) used by medin Medical Innovations GmbH see www.medin-medical.com/trademarks. © 2025 medin Medical Innovations GmbH. All rights reserved.

