


Nuflow

nasal cannula



Available sizes:

	Small	Medium	Large	X-Large
REF	1320	1330	1340	1350
Flow limits ¹	10,0 l/min	14 l/min	23 l/min	27 l/min
Outer diameter prong	2,3mm	2,7mm	3,2mm	3,7mm
Septum distance	3,7mm	4,3mm	4,8mm	5,3mm
Patient group ²	Newborns	Newborns	Infants	Children
Packing Unit	10 pieces	10 pieces	10 pieces	10 pieces
Fixation patches (Spare part)	1380 small	1380 small	1381 large	1381 large

Connectors

Tube connection: OD 10 mm or OD 15 mm available

Material

Cannula	Silicone, Zylar, TPE, PP
Fixation patch	Hydrocolloid
free of PVC, DEHP, natural rubber latex	

1 - See user manual for details
2 - Depends on size of patient nostril

Fixation patches (spare part)

REF 1380	Pair of Fixation-patches (small), (spare part for REF 1320/REF1330)
REF 1381	Pair of Fixation-patches (large), (spare part for REF 1340/1350)

Starter Sets

REF 1393	Nuflow starter set s/m 5x Nuflow nasal cannula small, 5x Nuflow nasal cannula medium
REF 1394	Nuflow starter set l/xl 5x Nuflow nasal cannula large 5x Nuflow nasal cannula x-large

Combinable devices

medinSINDI, medinCNOmini, medin-NC3, medinBLENDER



medin Medical Innovations GmbH
Adam-Geisler-Straße 1
82140 Olching, Germany
☎ +49 (0)8142 448 46 -0
✉ info@medin-medical.com
www.medin-medical.com

