





# REF 1207

Heated circuit set for Medijet® in combination with humidifier in use for incubator and baby warmer

clean - non steril



medinCNO®



medin-NC3®



Medijet®

## Specification

- Circuit gas supply inlet to humidifier (A):
  - Connectors both sides 22mm ID - female
  - Material tube: EVA, LDPE
  - Material connectors: K-Resin
- Circuit humidifier to Medijet (B):
  - Connector Medijet side: 7,4mm OD - male
  - Connector humidifier side: 22mm ID - female
  - Material tube: PP, EVA, LDPE
  - Material connectors: K-Resin, PP
- Pressure line (D):
  - Connectors both sides: Luer connectors
  - Material tube: PP
  - Material connectors: K-Resin
- Additional connectors:
  - TPR connector (C) to disconnect the front part of the circuit: 13mm ID - female (Medijet can be connected directly to the TPR connector)
  - Connector (E): 22mm OD - male; 15mm OD - male
  - Connector (F): 22mm OD - female; 22mm OD - male
  - Connector (G): 10mm OD - male; 10mm OD - male
  - Connector (H): 22mm OD - male; 22mm OD - male
  - Material connectors: E - H: K-Resin, C: TPR (Thermal plastic rubber)
- Heating:
  - The circuit humidifier to Medijet (B) includes a heating wire
  - Electric resistance 22,5 ohm
  - Electric connector: Cylindrical shape
  - Material: Brass, tin, PP
  - Two ports for the temperature sensors are included in the circuit
  - Material connectors: K-Resin, PP



medin Medical Innovations GmbH

Adam-Geisler-Straße 1

82140 Olching, Germany

+49 (0)8142 448 46 -0

info@medin-medical.com

[www.medin-medical.com](http://www.medin-medical.com)

