




1207 nCPAP circuit set


Technical specifications

For Medijet in combination with humidifier in use for incubator and baby warmer.

clean - non steril



medinCNO



medin-NC3



Medijet

Specifications

Circuit gas supply inlet to humidifier (A):

- Connectors both sides 22mm ID - female
- Material tube: EVA, LDPE
- Material connectors: K-Resin

Circuit humidifier to Medijet (B):

- Connector Medijet side: 7,4mm OD - male
- Connector humidifier side: 22mm ID - female
- Material tube: PP, EVA, LDPE
- Material connectors: K-Resin, PP

Pressure line (D):

- Connectors both sides: Luer connectors
- Material tube: PP
- Material connectors: K-Resin

Additional connectors:

- TPR connector (C) to disconnect the front part of the circuit: 13mm ID - female (Medijet can be connected directly to the TPR connector)
- Connector (E): 22mm OD - male; 15mm OD - male
- Connector (F): 22mm OD - female; 22mm OD - male
- Connector (G): 10mm OD - male; 10mm OD - male
- Connector (H): 22mm OD - male; 22mm OD - male
- Material connectors: E - H: K-Resin, C: TPR (Thermal plastic rubber)

Heating:

- The circuit humidifier to Medijet (B) includes a heating wire
- Electric resistance 22,5 ohm
- Electric connector: Cylindrical shape
- Material: Brass, tin, PP
- Two ports for the temperature sensors are included in the circuit
- Material connectors: K-Resin, PP



medin Medical Innovations GmbH

Adam-Geisler-Straße 1

82140 Olching, Germany

☎ +49 (0)8142 448 46 -0

info@medin-medical.com

www.medin-medical.com

REV07
19.05.2025

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. For all proprietary trademarks (®) and third-party trademarks (®) used by medin Medical Innovations GmbH see www.medin-medical.com/trademarks. © 2025 medin Medical Innovations GmbH. All rights reserved.