


# Nuflow<sup>®</sup>

## Nasal cannula set




- ✓ High degree of patient comfort
- ✓ Low noise level and system stagnation pressure
- ✓ Individual fixation, gentle to the skin
- ✓ Anatomically adapted shape of the prongs

Available sizes:

|                            | 1 Nasal cannula set <b>small</b> | 2 Nasal cannula set <b>medium</b> | 3 Nasal cannula set <b>large</b> | 4 Nasal cannula set <b>xlarge</b> |
|----------------------------|----------------------------------|-----------------------------------|----------------------------------|-----------------------------------|
| REF                        | 1320                             | 1330                              | 1340                             | 1350                              |
| Flow limits <sup>1</sup>   | 10,0 l/min                       | 14 l/min                          | 23 l/min                         | 27 l/min                          |
| Outer diameter prong       | 2,3 mm                           | 2,7 mm                            | 3,2 mm                           | 3,7 mm                            |
| Septum distance            | 3,7 mm                           | 4,3 mm                            | 4,8 mm                           | 5,3 mm                            |
| Patient group <sup>2</sup> | Premature                        | Premature, Newborns               | Infants                          | Children                          |
| Packing unit               | 10 pieces                        | 10 pieces                         | 10 pieces                        | 10 pieces                         |
| Fixation patches           | 1380 small                       | 1380 small                        | 1381 large                       | 1381 large                        |

1 - See user manual for details  
2 - Depends on size of patient nose



### Combinable devices

medinSINDI<sup>®</sup>, medinCNO<sup>®</sup>mini, medin-NC3<sup>®</sup>, medinBLENDER Ventilators with high-flow mode and overpressure protection. Other devices on request at the medin clinical application team.



Manufacturer:

medin Medical Innovations GmbH

Adam-Geisler-Straße 1

82140 Olching, Germany

☎ +49 (0)8142 448 46 -0

info@medin-medical.com

**[www.medin-medical.com](http://www.medin-medical.com)**

REV02  
16.05.2019

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. For all proprietary trademarks and third-party trademarks used by medin Medical Innovations GmbH see [www.medin-medical.com/trademarks](http://www.medin-medical.com/trademarks).  
© 2018 medin Medical Innovations GmbH. All rights reserved.

Nuflow

